

PVMS OF MEDICAL EQUIPMENT

Clinical Specialty	Radiological Equipment
Generic Name	X-RAY FILM PROCESSOR
Clinical Purpose	It helps in automatic processing of x-ray films providing better quality images and consuming less time.

TECHNICAL SPECIFICATIONS

High throughput free Standing type
Continuous roller transport
Sheet film, 10x10cm-35x43cm sizes (14"x17")
High through put more than 110 sheets (35 x 43cm)
Forced hot air circulation method
Piping installed in the DEV and FIX racks.
Processing solution temperature; Controlled by the temperature control tank, with the thermistor monitoring, with the heater heating and with the wash water cooling. Drying temperature; Controlled automatically according to a fixed temperature setting.

APPROVED PVMS